

OSCE Code 7457-350 Theatre Practice 04 (Pack drape) for Unit 310

This station covers the following RCVS Day One Skills: 9.2

This OSCE will be used to assess the awards indicated

Award	Award Reference	Pathway
7457 – Level 3 Diploma in Veterinary Nursing	600/6052/6	Small Animal
7457 – Level 3 Diploma in Veterinary Nursing	600/6052/9	Equine

Scenario

Prior to a surgical procedure, you have been requested to:

- a. Select equipment and place on the correct section of the grid provided
- b. Pack and label the drape indicated by the examiner, ready for sterilisation using an autoclave today

Note: The equipment has been previously cleaned

Methodology: You will be expected to:

1.	From the equipment provided correctly fill 9 spaces on the grid
2.	Fold the drape using a recognised method
3.	Select appropriate sized autoclave packaging material
4.	Select TST indicator strip
5.	Insert a TST inside the folded drape
6.	Insert drape into autoclave packaging
7.	Remove excess air from the bag
8.	Secure the package
9.	Double wrap in appropriate sized autoclave packaging material
10.	Remove excess air from bag
11.	Secure the package
12.	Select permanent marker pen
13.	Label package with name of contents
14.	Label package with quantity of drapes
15.	Label package with own initials
16.	Label package with date of autoclaving or date for re-autoclaving
17.	Write on appropriate portion of packaging without piercing bag (ie flap, plastic or Bowie- Dick)
18.	Drape packed appropriately ready for sterilisation with the correct inclusion of the TST strip