

OSCE Code 7457-350 Theatre Practice 03 (Pack instrument) for Unit 310

This station covers the following RCVS Day One Skills: 9.2

This OSCE will be used to assess the awards indicated

| Award | Award Reference | Pathway |
|--|-----------------|--------------|
| 7457 – Level 3 Diploma in Veterinary Nursing | 600/6052/6 | Small Animal |
| 7457 – Level 3 Diploma in Veterinary Nursing | 600/6052/9 | Equine |

Scenario

Prior to a surgical procedure, you have been requested to identify and prepare equipment for sterilisation

- a. Select the required items and place them on the correct section of the grid provided
- b. Pack and label the equipment indicated by the examiner, ready for sterilisation using an autoclave today

Note: The equipment has been previously cleaned

| Met | Methodology: You will be expected to: | | |
|-----|--|--|--|
| 1. | From the equipment provided correctly fill 9 spaces on the grid | | |
| 2. | Cover the sharp end of the equipment securely | | |
| 3. | Select appropriate sized autoclave packaging material | | |
| 4. | Insert the equipment into packaging | | |
| 5. | Equipment to be inserted the correct way round (ie handles are presented first when packaging is opened) | | |
| 6. | Remove excess air from the bag | | |
| 7. | Secure the package | | |
| 8. | Double wrap in appropriate sized autoclave packaging material | | |
| 9. | Remove excess air from the bag | | |
| 10. | Secure the package | | |
| 11. | Select permanent marker pen | | |
| 12. | Label package with correct equipment name | | |
| 13. | Label package with own initials | | |
| 14. | Label package with date of autoclaving or date for re-autoclaving | | |
| 15. | Write on appropriate portion of packaging without piercing bag (ie flap, plastic or Bowie-Dick) | | |
| 16. | Item of equipment packed appropriately ready for sterilisation | | |